

Midmark – Mobile Dental Delivery System

Conçu pour une plus grande efficacité dans votre pratique vétérinaire

Les soins de santé bucco-dentaire sont importants pour vos patients et votre cabinet. Le système dentaire mobile Midmark de nouvelle génération contribue à améliorer votre façon de travailler en rationalisant le flux de travail dentaire grâce à des outils et des ressources intégrés qui fournissent tout ce dont vous avez besoin pour effectuer vos traitements dentaires.

- Conception compacte
- Rangement pratique
- Grande surface de travail
- Poignées intégrées
- Commandes conviviales
- Gestion du cordon
- Tubes de pièces à main lisses
- Indicateur LED de niveau d'eau
- Système de bouteille d'eau à connexion/déconnexion rapide
- Points d'accès faciles pour l'entretien
- Liste de contrôle de l'entretien à bord
- Pièce à main et détartréur à LED haute vitesse de série
- Deux pièces à main haute vitesse en option
- Aspiration intégrée en option



Caractéristiques techniques

Weight & Dimensions

Overall depth: 18.6" (47.24 cm)

Overall width: 17.4" (44.20 cm)

Overall height: 33.8" (85.85 cm)

Weight: 130 lbs

Compresso:

Oil-free 1/2 hp motor

Working pressure: 70-100 p.s.i.

Flow rate: 2.2 CFM at 70 p.s.i.

Air reservoir: 1 gallon

Manual air tank purge

Water Module

2 - 1 liter reservoirs (Toggled)

Pressurized system: 30-40 psi

Manual handpiece flush system

External quick connects for water and air

Electrical Ratings

Voltage:

- 115 VAC, (models: -001, -002, -003, -005)
- 230 VAC (model: -004)

Frequency: 50/60 Hz

Circuit breaker protection:

- 6A, (models: -001, -002, -003, -005)
- 3A (model: -004)

Powercord length: 10'

Model Numbers

AH9000-001 Midmark Mobile Dental Delivery System, 4-Position Integrated LED High-Speed Handpiece and Scaler

AH9000-002 Midmark Mobile Dental Delivery System, Swivel, 4 Position Integrated LED High-Speed Handpiece and Scaler

AH9000-003 Midmark Mobile Dental Delivery System, Swivel, 5 Position with 2 Integrated LED High-Speed Handpieces and a Scaler

AH9000-004 Midmark Mobile Dental Delivery System, Integrated 4-Position Export

AH9000-005 Midmark Mobile Dental Delivery System, 4-Position Integrated LED High-Speed Handpiece and Scaler, CANADA

