

Midmark – Mobile Dental Delivery System

Ontworpen voor meer efficiëntie in uw dierenartsenpraktijk

Mondverzorging is belangrijk voor uw patiënten en uw praktijk. Het mobiele tandheelkundige systeem van Midmark van de volgende generatie helpt de manier waarop u werkt te verbeteren door de tandheelkundige workflow te stroomlijnen met geïntegreerde tools en hulpmiddelen die alles bieden wat u nodig hebt om tandheelkundige behandelingen uit te voeren.

- Compact ontwerp
- Handige opslag
- Groot werkoppervlak
- Geïntegreerde handgrepen
- Gebruiksvriendelijke bediening
- Snoerbeheer
- Gladde handstukbuizen
- LED waterpeilindicator
- Snel aan-/afkoppelbaar waterflessysteem
- Gemakkelijke toegangspunten voor onderhoud
- Controlelijst voor onderhoud aan boord
- Standaard LED-handstuk met hoge snelheid en LED-scaler
- Optioneel dubbele hogesnelheidshandstukken
- Optionele geïntegreerde afzuiging



Technische kenmerken

Weight & Dimensions

Overall depth: 18.6" (47.24 cm)

Overall width: 17.4" (44.20 cm)

Overall height: 33.8" (85.85 cm)

Weight: 130 lbs

Compresso:

Oil-free 1/2 hp motor

Working pressure: 70-100 p.s.i.

Flow rate: 2.2 CFM at 70 p.s.i.

Air reservoir: 1 gallon

Manual air tank purge

Water Module

2 - 1 liter reservoirs (Toggled)

Pressurized system: 30-40 psi

Manual handpiece flush system

External quick connects for water and air

Electrical Ratings

Voltage:

- 115 VAC, (models: -001, -002, -003, -005)
- 230 VAC (model: -004)

Frequency: 50/60 Hz

Circuit breaker protection:

- 6A, (models: -001, -002, -003, -005)
- 3A (model: -004)

Powercord length: 10'

Model Numbers

AH9000-001 Midmark Mobile Dental Delivery System, 4-Position Integrated LED High-Speed Handpiece and Scaler

AH9000-002 Midmark Mobile Dental Delivery System, Swivel, 4 Position Integrated LED High-Speed Handpiece and Scaler

AH9000-003 Midmark Mobile Dental Delivery System, Swivel, 5 Position with 2 Integrated LED High-Speed Handpieces and a Scaler

AH9000-004 Midmark Mobile Dental Delivery System, Integrated 4-Position Export

AH9000-005 Midmark Mobile Dental Delivery System, 4-Position Integrated LED High-Speed Handpiece and Scaler, CANADA

